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
Exempel:

	Vänster hand					Höger hand			
Former									
	A	B	C			A	B	C	
Ø15	✓	✓	✓	1		✓	✓	x	0
Ø8	✓	✓	✓	1					
Ø5	✓	✓	✓	1					
Ytor									
	D	E	F			D	E	F	
15 mm	✓	✓	✓	1		✓	✓	✓	1
8 mm	✓	✓	✓	1		✓	x		0
4 mm	✓	✓	✓	1					
	Summa			6		Summa			1

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Kommentarer:

STI Literature:

In chronological order, includes a selection of studies that have utilized the STI as an outcome measure

Rosén B. Recovery of sensory and motor function after nerve repair. A rationale for evaluation. *J Hand Ther* 1996;9:315–327.

Rosén B, Lundborg G. A new tactile gnosis instrument in sensibility testing. *J Hand Therapy* 1998;11:251-257.

Rosén B, Lundborg G. A model instrument for the documentation of outcome after nerve repair. *J Hand Surg*, 2000;25A:535-543.

Rosén B, Jerosch-Herold C. Comparing the responsiveness over time of two tactile gnosis tests: two-point discrimination and the STI-test. *Br J Hand Ther*, 2000;5(4):114-119.

Rosén B, Lundborg G. The long term recovery curve in adults after median or ulnar nerve repair: a reference interval. *Journal of Hand Surgery (British and European Volume)*. 2001; 26(3):196-200.

Rosén B. Inter-tester reliability of a tactile gnosis test: the STI-test., *British J Hand Ther* 2003;8(3):98-101

Rosén B, Lundborg G. A new model instrument for outcome after nerve repair. *Hand clinics*. 2003; 19(3):463-70.

Jerosch-Herold C. Assessment of sensibility after nerve injury and repair: a systematic review of evidence for validity, reliability and responsiveness of tests. *The Journal of Hand Surgery: British & European Volume*. 2005 Jun 30;30(3):252-64.

Svensson E, Häger-Ross C. Hand function in Charcot Marie Tooth : test retest reliability of some measurements. *Clinical Rehabilitation*. 2006;20(10)

Vordemvenne T, Langer M, Ochman S, Raschke M, Schult M. Long-term results after primary microsurgical repair of ulnar and median nerve injuries: A comparison of common score systems. *Clinical neurology and neurosurgery*. 2007; 109(3):263-71.

Eklund E, Svensson E, Häger-Ross C. Hand function and disability of the arm, shoulder and hand in Charcot-Marie-Tooth disease. *Disabil Rehabil*. 2009;31(23):1955-62.

Jerosch-Herold C, Shepstone L, Miller L, Chapman P. The responsiveness of sensibility and strength tests in patients undergoing carpal tunnel decompression 2011, 12:244

Post R, de Boer KS, Malessy MJ. Outcome following nerve repair of high isolated clean sharp injuries of the ulnar nerve. *PLoS One*. 2012 Oct 17;7(10):e47928.

Miller LK, Chester R, Jerosch-Herold C. Effects of sensory reeducation programs on functional hand sensibility after median and ulnar repair: a systematic review. *Journal of Hand Therapy*. 2012 Sep 30;25(3):297-307.

Post R, de Boer KS, Malessy MJA (2012) Outcome following Nerve Repair of High Isolated Clean Sharp Injuries of the Ulnar Nerve. PLOS ONE 7(10): e47928. <https://doi.org/10.1371/journal.pone.0047928>

Rosén B, Björkman A, Boeckstyns M. Differential recovery of touch thresholds and discriminative touch following nerve repair with focus on time dynamics. *Hand Therapy*. 2014;1758998314538004.

van Gils W, Reinders-Messelink HA, Smit-Klaaij F, Bongers RM, Dijkstra PU, van der Sluis CK. Sensibility of the stump in adults with an acquired major upper extremity amputation. *Archives of physical medicine and rehabilitation*. 2013;94(11):2179-85.

Rosén B, Vikström P, Turner S, McGrouther DA, Selles RW, Schreuders TA, Björkman A. Enhanced early sensory outcome after nerve repair as a result of immediate post-operative re-learning: a randomized controlled trial. *Journal of Hand Surgery (European Volume)*. 2015;40(6):598-606.

Ekstrand E, Lexell J, Brogårdh C. Test-retest reliability of the Shape/Texture Identification testTM in people with chronic stroke. *Clin Rehabil*. 2016 Nov;30(11):1120-1127.

Jerosch-Herold, C., Houghton, J., Miller, L., & Shepstone, L. (2016). Does sensory relearning improve tactile function after carpal tunnel decompression? A pragmatic, assessor-blinded, randomized clinical trial. *The Journal of Hand Surgery, European Volume*, 41(9), 948–956.

If you have published a study utilizing the STI, please let us know, we may include it in future literature lists.

STI₂ Instructions