Bruksanvisning – Catell Armsling Double

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230416:4

Instruction for use – Catell Armsling Double

Catell recommend that the product is fitted and adapted by a licensed or otherwise qualified health professionals. Users with diabetes, sensitive skin, known contact allergies, impaired sensibility or poor circulation must be extra observant and check for pinching, chafing or edema when using the brace. Always check material content below. Any serious incidents occurring relates to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident. Contact your prescriber or Catell if questions about the product occur.

Area of use

The brace supports the users arm and shoulder after injury or surgery.

Application

1. First adjust the arm pouch with its wide hook closure so that it fits the intended arm.
2. Open the hook and loop closures and place the forearm in the pouch. Picture 1
3. Pull the long straps over the shoulders, cross them behind your back, place around the forearm and attach the hook closure on the arm pouch. Thoroughly adjust the correct strap length for optimal support. Picture 2a, 2b
4. In need of additional fixation of the arm to the body, use the long wide additional strap. Picture 3
5. Make sure that all straps are properly secured.

When necessary adjust the brace. It is important that the brace fits comfortably without pinching, chafing or affecting blood circulation.

Material

Nylon, polyester. Free from latex.

Washing instructions

Close the hook and loop straps. Follow the washing instructions below. Use a washing pouch. Let the brace drip dry.

Size

Size	*Dimensions	Size	*Dimensions	Size	*Dimensions
S	28 - 33 cm	M	33 - 38 cm	L	38 - 43 cm

*Measure length of forearm. Measurements should be seen as guidelines. The product is bilateral. Single patient use.

Recycling

Separate plastic and textile. Sort according to local regulations.

This device conforms to the requirements of the Medical Device Regulation (EU) 2017/745 product class I.

Manufacturer Catell AB
Örkroken 23, SE-138 40 Älta, Sweden | Phone: +46 8 727 95 30 | info@catell.se | www.catell.se

230416:4